

Restore True Color, Enjoy Color Matching



YS3060 Accuracy Spectrophotometer

YS3060 is independently developed by 3nh, who has completed intellectual property rights. With full set of light sources, dual apertures (8&4mm), SCI/SCE, UV, USB/Bluetooth, it has high accuracy and large storage, very suitable for lab color analysis and transmission. Accuracy, stability, powerful functions make it in the leading position of color industry.







USB/Bluetooth



LED/UV light



Camera Locating

















Product Features -----

- 1.D/8 geometrical optics, conforms with CIE No.15,GB/T 3978,GB2893 GB/T 18833, ISO7724/1, ASTM E1164, DIN5033 Teil7
- 2. Combined long life LED light source, including/excluding UV
- 3. Switchable Φ8mm and Φ4mm 2 apertures, adapt to more samples
- 4.USB/Bluetooth 4.0, dual modes, widely useful
- 5. Support both SCI and SCE modes measurement
- 6. Camera locating function, better position
- 7. Super stain-resistant and stable white calibration plate
- 8. Large capacity storage space, over 20,000 measurement data
- 9. With full illuminants and many color indexes for color measurement
- 10.PC software has a powerful function extension



Application Industries --















Automobile

Leather

Plastics

Paint

Food stuff

Laboratory

Others

Specification Parameters

Model: YS3060 Spectrophotometer

Optical Geometry: Reflect: di:8°, de:8°(diffused illumination, 8-degree viewing angle)

Integrating Sphere Size: 48mm

Light Source: Combined LED Light, UV Light **Spectrophotometric Mode:** Concave Grating

Sensor: 256 Image Element Double Array CMOS Image Sensor

Wavelength Range: 400-700nm Wavelength Interval: 10nm Semiband Width: 10nm

Measured Reflectance Range: 0-200%

Measuring Aperture: Dual Aperture: 10mm/8mm & 5mm/4mm

Specular Component: SCI&SCE

Color Space: CIE Lab, XYZ, Yxy, LCh, CIE LUV, Hunter LAB

 $\textbf{Color Difference Formula:} \ \Delta \text{E*ab}, \ \Delta \text{E*uv}, \ \Delta \text{E*94}, \ \Delta \text{E*cmc}(2:1), \ \Delta \text{E*cmc}(1:1), \ \Delta \text{E*00v}, \ \Delta \text{E}(\text{Hunter})$

Other Colorimetric Index: WI(ASTM E313, CIE/ISO, AATCC, Hunter),

YI(ASTM D1925, ASTM 313), TI(ASTM E313, CIE/ISO),

 $Metamerism\ Index\ MI,\ Staining\ Fastness,\ Color\ Fastness,$

Color Strength, Opacity, 8° Glossiness

Illuminant: D65, A, C, D50, D55, D75, F1, F2, F3, F4, F5, F6, F7, F8, F9, F10, F11, F12 Displayed Data: Spectrogram/Values, Samples Chromaticity Values, Color Difference

Values/Graph, PASS/FAIL Result, Color Offset

Observer Angle:2°/10° Measuring Time: 2.6s

Repeatability:MAV/SCI: ∆E*≤0.03

Inter-instrument Error: MAV/SCI: ∆E*≤0.15

 $\textbf{Measurement Mode:} \ \mathsf{Single Measurement}, \mathsf{Average Measurement}$

Locating Method: Camera Locating

Battery: Li-ion battery. 5000 measurements within 8 hours

Dimension: L*W*H=184*77*105mm

Weight: 600g

 $\textbf{Illuminant Life Span:}\ 5\ \text{years, more than}\ 3\ \text{million times measurements}$

Display: 3.5-inch TFT color LCD, Capacitive Touch Screen

Data Port: USB, Bluetooth 4.0

Data Storage: Standard 1000 Pcs, Sample 28000 Pcs

Language: English, Chinese

 $\textbf{Operating Environment:} \ 0\sim40^{\circ}\text{C},\ 0\sim85\%\text{RH (no condensing)},\ Altitude < 2000\text{m}$

Storage Environment: -20~50°C, 0~85%RH (no condensing)

Standard Accessory: Power Adapter, Built-In Li-ion Battery, User Guide, PC

Software, White and Black Calibration Cavity, Dust Cover

Optional Accessory: Micro Printer, Powder Test Box



SHENZHEN ThreeNH TECHNOLOGY CO., LTD.

Address: F/6, Block 5B, Skyworth Inno Valley, Tangtou 1st Road, Shiyan, Baoan District, Shenzhen, P.R. China

Service Hotline: 400-666-2522 Email: 3nh@3nh.com
Tel: 86-0755-26508999 Fax: 86-0755-27190609